


Summary

Name	UBP34_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q70CQ2
Host	E.coli
Species	Human
Molecular Mass	39 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EQLINRTKHV QQLSDTEES MQGSSDETAN SGEDGSSGPG SSSGHS DGSS NEVNSSHASQ SAGSPGSEVQ SEDIADIEAL KEDEDDDHG HNPPKSSCGT DLRNRKLESQ AGICLGDSQG MSERNGTSSG TGKDLVFNTE SLPSVDNRMR MLDACSHSED PEHDISGEMN ATHIAQGSQE SCITRTGDFL GETIGNELFN CRQFIGPQHH HHHHHHHHHH DGHMVDDMLS ADDVSCSSSQ VSAKSEKNMA DFDGEESGCE EELVQINSHA ELTSHLQQHL PNLASIYHEH LSQGPVVHKH QFNSNAVTDI NLDNVCKKGN TLLWDIVQDE DAVNLSEGLI

image



Background

Alternative Names	Ubiquitin carboxyl-terminal hydrolase 34 (EC 3.4.19.12) (Deubiquitinating enzyme 34) (Ubiquitin thioesterase 34) (Ubiquitin-specific-processing
--------------------------	---

Product Name: Recombinant Human USP34 Protein
Catalog #: PEH9234



protease 34)

Note

For Research Use Only , Not for Diagnostic Use