


Product Name: Recombinant Human RAPH1 Protein
Catalog #: PEH9235



Summary

Name	RAPH1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q70E73
Host	E.coli
Species	Human
Molecular Mass	41 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	HGAEEDSDKE DQDLDKMFGA WLGELDKLTQ SLDSDKPMEP VKRSPLRQET NMANFSYRFS IYNLNEALNQ GETVDLDALM ADLCSIEQEL SSIGSGNSKR QITETKATQK LPVSRHTLKH GTLKGLSSSS NRIAKPSHAS YSLDDVTAQL EQASLSMDEA AQQSVLEDTK PLVTNQHRRT ASAGTVSDAE VHSISNSSHS SITSAASSMD SLDIDKVTRP QELDLTHQGQ PITEEEQAAK LKAEKIRVAL EKIKEAQVKK LVIRVHMSDD SSKTMMVDER QTVRQVLDNL MDKSHCGYSL DWSLVETVSE LQMERIFEDH ENLVENLLNW TRDSQNKLIF

image



Background

Alternative Names Ras-associated and pleckstrin homology domains-containing protein 1 (RAPH1) (Amyotrophic lateral sclerosis 2 chromosomal region candidate)

Product Name: Recombinant Human RAPH1 Protein
Catalog #: PEH9235



gene 18 protein) (Amyotrophic lateral sclerosis 2 chromosomal region
candidate gene 9 protein) (Lamellipodin) (Proline-rich EVH1 ligand 2)
(PREL-2) (Protein RMO1)

Note

For Research Use Only , Not for Diagnostic Use