


Product Name: Recombinant Human MICAL3 Protein
Catalog #: PEH9241



Summary

Name	MICA3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q7RTP6
Host	E.coli
Species	Human
Molecular Mass	29 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	DSGKHRVLPV VRAELQLRRQ LSFSESDLS SDDVLEKSSQ KSRREPRTYT EEELNAKLTR RVQKAARRQA KQEELKRLHR AQIIQRQLQQ VEERQRRLEE RGVAVEKALR GEAGMGKKDD PKLMQEWFKL VQEKNAMEVRY ESELMIFARE LELEDQRSRL QQELRERMAV EDHLKTEEEL SEEKQILNEM LEVVEQRDSL VALLEEQLRL EREEDKDLEA AMLSKGFSLN

image



Background

Alternative Names	[F-actin]-monooxygenase MICAL3 (EC 1.14.13.225) (Molecule interacting with CasL protein 3) (MICAL-3)
--------------------------	--

Note

Product Name: Recombinant Human MICAL3 Protein
Catalog #: PEH9241



For Research Use Only , Not for Diagnostic Use