


Product Name: Recombinant Human TTC21B Protein
Catalog #: PEH9244



Summary

Name	TT21B_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q7Z4L5
Host	E.coli
Species	Human
Molecular Mass	12 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	DNETVGGEVF ENLDGDLGNS TEKQESVQLA VRTAEKLLKE LKPQTVQGHV QLRIMENYCL MATKQKSNVE QALNTFTEIA

image



Background

Alternative Names Tetratricopeptide repeat protein 21B (TPR repeat protein 21B) (Intraflagellar transport 139 homolog)

Note

For Research Use Only , Not for Diagnostic Use