


Product Name: Recombinant Human SRGAP1 Protein
Catalog #: PEH9246



Summary

Name	SRGP1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q7Z6B7
Host	E.coli
Species	Human
Molecular Mass	41 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EHGTL EEVDQ DAGTEPHTSE DECEPIEAIA KFDYVGRSAR ELSFKKGASL LLYHRASEDW WEGRHNGIDG LVPHQYIVVQ DMDDTFSDTL SQKADSEASS GPVTEDKSSS KDMNSPTDRH PDGYLARQRK RGEPPPPVRR PGRTSDGHCP LHPPHALSNS SVDLGSPSLA SHPRGLLQNR GLNNDSPERR RRPGHGSLTN ISRHDSLKKI DSPPIRRSTS SGQYTGFNHDH KPLDPETIAQ DIEETMNTAL NELRELERQS TAKHAPDVVL DTLEQVKNSP TPATSTESLS PLHNVALRSS EPQIRRSTSS SSDTMSTFKP MVAPRMGVQL KPPALRPKPA VLPKTNPTIG PAPPPQGPTD

image



Background

Alternative Names SLIT-ROBO Rho GTPase-activating protein 1 (srGAP1) (Rho GTPase-activating

Product Name: Recombinant Human SRGAP1 Protein
Catalog #: PEH9246



protein 13)

Note

For Research Use Only , Not for Diagnostic Use