


Product Name: Recombinant Human NECAB2 Protein
Catalog #: PEH9247



Summary

Name	NECA2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q7Z6G3
Host	E.coli
Species	Human
Molecular Mass	18 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GTAVILDIFR RADKNDDGKL SLEEFQLFFA DGVLNKELE DLFHTIDSDN TNHVDTKELC DYFVDHMGDY EDVLASLETL NHSVLKAMGY TKKVYEGGSN VDQFVTRFLL KETANQIQSL LSSVESAVEA IEEQTSQLRQ

image



Background

Alternative Names	N-terminal EF-hand calcium-binding protein 2 (EF-hand calcium-binding protein 2) (Neuronal calcium-binding protein 2) (Synaptotagmin-interacting protein 2) (Stip-2)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human NECAB2 Protein
Catalog #: PEH9247

