


Product Name: Recombinant Human ADAT2 Protein
Catalog #: PEH9249



Summary

Name	ADAT2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q7Z6V5
Host	E.coli
Species	Human
Molecular Mass	19 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	ASGACSVSAE ETEKWMEEAM HMAKEALENT EVPVGCLMVY NNEVVGKGRN EVNQTKNATR HAEMVAIDQV LDWCRQSGKS PSEVFEHTVL YVTVEPCIMC AAALRLMKIP LVVYGQNER FGGCGSVLNI ASADLPNTGR PFQCIPGYRA

image



Background

Alternative Names	tRNA-specific adenosine deaminase 2 (EC 3.5.4.33) (Deaminase domain-containing protein 1) (tRNA-specific adenosine-34 deaminase subunit ADAT2)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human ADAT2 Protein
Catalog #: PEH9249

