


Summary

Name	UBP37_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q86T82
Host	E.coli
Species	Human
Molecular Mass	20 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KETSRQLSYS DNQASAKRGS LETKDDIPFR KVLGNPGRGS IKTVAGSGIA RTIPSLTSTS TPLRSGLEN RTEKRKRMIS TGSELNEDYP KENDSSSNK AMTDPSRKYL TSSREKQLSL KQSEENRTSG LLPLQSSSFY GSRAGSKEHS SGGTNLDRTN VSSQTPSAKR

image



Background

Alternative Names	Ubiquitin carboxyl-terminal hydrolase 37 (EC 3.4.19.12) (Deubiquitinating enzyme 37) (Ubiquitin thioesterase 37) (Ubiquitin-specific-processing protease 37)
--------------------------	--

Note

Product Name: Recombinant Human USP37 Protein
Catalog #: PEH9250



For Research Use Only , Not for Diagnostic Use