


Summary

Name	PB1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q86U86
Host	E.coli
Species	Human
Molecular Mass	39 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs.Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KNTDSEDSR AEDNFNLEKE KEDVPVEMSN GEPGCHYFEQ LHYNDMWLKV GDCVFIKSHG LVRPRVGRIE KVWVRDGAAY FYGPIFIHPE ETEHEPTKMF YKKEVFLSNL EETCPMTCIL GKCAVLSFKD FLSCRPTTEIP ENDILLCESR YNESDKQMKK FKGLKRFSLS AKVVDDEIYY FRKPIVPQKE PSQLLEKKIQ LLEAKFAELE GGDDDIEMG EEDSEVIEPP SLPQLQTPLA SELDLMPYTP PQSTPKSAKG SAKKEGSKRK INMSGYILFS SEMRAVIKAQ HPDYSFGELS RLVGTEWRNL ETAKKAEYEE

image



SDS



WB

Background

Alternative Names	Protein polybromo-1 (hPB1) (BRG1-associated factor 180) (BAF180) (Polybromo-1D)
--------------------------	---

Product Name: Recombinant Human PBRM1 Protein
Catalog #: PEH9251



Note

For Research Use Only , Not for Diagnostic Use