


Product Name: Recombinant Human NAA40 Protein
Catalog #: PEH9254



Summary

Name	NAA40_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q86UY6
Host	E.coli
Species	Human
Molecular Mass	28 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KKQKRLEERA AMDAVCAKVD AANRLGDPLE AFPVFKKYDR NGLNVSIECK RVSGLEPATV DWAFDLTKTN MQTMYEQSEW GWKDREKREE MTDDRAWYLI AWENSSVPVA FSHFRFDVEC GDEVLYCYEV QLESKVRRKG LGKFLIQILQ LMANSTQMKK VMLTVFKHNH GAYQFFREAL QFEIDDSSPS MSGCCGEDCS YEILSRRTKF GDSHHSHAGG

image



Background

Alternative Names	N-alpha-acetyltransferase 40 (EC 2.3.1.257) (N-acetyltransferase 11) (N-alpha-acetyltransferase D) (NatD) (hNatD) (Protein acetyltransferase 1)
--------------------------	---

Note

Product Name: Recombinant Human NAA40 Protein
Catalog #: PEH9254



For Research Use Only , Not for Diagnostic Use