


Product Name: Recombinant Human DDX42 Protein
Catalog #: PEH9259



Summary

Name	DDX42_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q86XP3
Host	E.coli
Species	Human
Molecular Mass	17 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SAGAQGVNNT ASGNSREGT GGSNGKRERY TENRGSSRHS HGETGNRHSD SPRHGDGGRH GDGYRHPESS SRHTDGHRHG ENRHGGSAGR HGENRGANDG RNGESRKEAF NRESKMEPKM EPKVDSSKMD KVDSKTDKTA

image



Background

Alternative Names	ATP-dependent RNA helicase DDX42 (EC 3.6.4.13) (DEAD box protein 42) (RNA helicase-like protein) (RHELP) (RNA helicase-related protein) (RNAHP) (SF3b DEAD box protein) (Splicing factor 3B-associated 125 kDa protein) (SF3b125)
--------------------------	---

Note

Product Name: Recombinant Human DDX42 Protein
Catalog #: PEH9259



For Research Use Only , Not for Diagnostic Use