


Product Name: Recombinant Human TRIM50 protein
Catalog #: PEH9260



Summary

Name	TRIS0_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q86XT4
Host	E.coli
Species	Human
Molecular Mass	13 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	RMKEELAALI SELKQEQKKV DELIAKLVNN RTRIVNESDV FSWVIRREFQ ELHHLVDEEK ARCLEGIGGH TRGLVASLDM QLEQAQGTRE RLAQAECVLE

image



Background

Alternative Names E3 ubiquitin-protein ligase TRIM50 (EC 2.3.2.27) (RING-type E3 ubiquitin transferase TRIM50) (Tripartite motif-containing protein 50)

Note

For Research Use Only , Not for Diagnostic Use