


Product Name: Recombinant Human MSRB3 Protein
Catalog #: PEH9265



Summary

Name	MSRB3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8IXL7
Host	E.coli
Species	Human
Molecular Mass	18 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	CKVVFSSQQEL RKRLTPLQYH VTQEKGTESA FEGEYTHHKD PGIYKCVVCG TPLFKSETKF DSGSGWPSFH DVINSEAITF TDDFSYGMHR VETSCSQCGA HLGHIFDDGP RPTGKRYCIN SAALSFTPAD SSGTAEGGSG VASPAQADKA

image



Background

Alternative Names Methionine-R-sulfoxide reductase B3 (MsrB3) (EC 1.8.4.12) (EC 1.8.4.14)

Note

For Research Use Only , Not for Diagnostic Use