


Product Name: Recombinant Human CCDC65 Protein
Catalog #: PEH9266



Summary

Name	DRC2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8IXS2
Host	E.coli
Species	Human
Molecular Mass	26 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SLAKDLSEAE EQYAHALRSH LHNVDQLLAL QRHRLSLLEE SYNMELEALT KEFETERKTI IDQHEKEIHY LQDIFMAMEQ NYIDSEYESK LEFQSMWNDL KNMNLEEKHF LRLHLENRVE DLWRKFQDVL KNYTDATEDR KAAFETLQVK DEKSSKEIEV QMKKIQLQD AITISKGKIM IHSRESEDEN

image



Background

Alternative Names	Dynein regulatory complex subunit 2 (Coiled-coil domain-containing protein 65) (Testis development protein NYD-SP28)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human CCDC65 Protein
Catalog #: PEH9266

