


**Product Name: Recombinant Human RPAP2 Protein**  
**Catalog #: PEH9267**



## Summary

|                              |                                                                                                                                                                                                                                                                                                                                                                                                                                 |
|------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Name</b>                  | RPAP2_HUMAN                                                                                                                                                                                                                                                                                                                                                                                                                     |
| <b>Purity</b>                | Greater than 85% as determined by reducing SDS-PAGE                                                                                                                                                                                                                                                                                                                                                                             |
| <b>Endotoxin level</b>       | Please contact with the lab for this information                                                                                                                                                                                                                                                                                                                                                                                |
| <b>Tag</b>                   | C-His                                                                                                                                                                                                                                                                                                                                                                                                                           |
| <b>Accession #</b>           | Q8IXW5                                                                                                                                                                                                                                                                                                                                                                                                                          |
| <b>Host</b>                  | E.coli                                                                                                                                                                                                                                                                                                                                                                                                                          |
| <b>Species</b>               | Human                                                                                                                                                                                                                                                                                                                                                                                                                           |
| <b>Molecular Mass</b>        | 43 kDa                                                                                                                                                                                                                                                                                                                                                                                                                          |
| <b>Buffer</b>                | Phosphate buffered saline                                                                                                                                                                                                                                                                                                                                                                                                       |
| <b>Shipping</b>              | The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.                                                                                                                                                                                                                                                                                                              |
| <b>Stability&amp;Storage</b> | Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.                                                                                                                                                                                                                  |
| <b>Sequence</b>              | IPKTPVWVRE EERHPDFQLL KEEQSGHSGE EVQLCSKAIK TSDIDNPSHF<br>EKQYESSSSS THSDSSSDNE QDFVSSILPG NRPNSTNIRP QLHQKSIMKK<br>KAGHKANSKH KDKEQTVVDV TEQLGDCKLD SQEKDATCEL PLQKVNTQSS<br>SNSTLPERLK ASENSESEYS RSEITLVGIS KKSAEHFKRK FAKSNQVSRS<br>VSSSVQVCPE VGKRNLLKVL KETLIEWKTE ETLRFLYGQN YASVCLKPEA<br>SLVKEELDED DIISDPDSHF PAWRESQNSL DESLPFRGSG TAIKPLPSYE<br>NLKKETEKLN LRIREFYRGR YVLGEETTKS QDSEHDSTF PLIDSSSQNQ<br>IRKRIVLEKL |

## image



## Background

|                          |                                                                           |
|--------------------------|---------------------------------------------------------------------------|
| <b>Alternative Names</b> | Putative RNA polymerase II subunit B1 CTD phosphatase RPAP2 (EC 3.1.3.16) |
|--------------------------|---------------------------------------------------------------------------|

**Product Name: Recombinant Human RPAP2 Protein**  
**Catalog #: PEH9267**



---

(RNA polymerase II-associated protein 2)

**Note**

For Research Use Only , Not for Diagnostic Use