


Product Name: Recombinant Human SMAP1 Protein
Catalog #: PEH9270



Summary

Name	SMAP1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8IYB5
Host	E.coli
Species	Human
Molecular Mass	29 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	QKLNEQHQLI LSKLLREEDN KYCADCEAKG PRWASWNIGV FICIRCAGIH RNLGVHISRV KSVNLDQWTA EQIQCMQDMG NTKARLLYEA NLPENFRRPQ TDQAVEFFIR DKYEKKKYD KNAIAITNIS SSDAPLQPLV SSPSLQAAVD KNKLEKEKEK KKEEKKREKE PEKPAKPLTA EKLQKKDQQL EPKKSTSPKK AAEPTVDLLG LDGPAVAPVT NGNTTVPLN DDLDFGPMI

image



Background

Alternative Names Stromal membrane-associated protein 1

Note

Product Name: Recombinant Human SMAP1 Protein
Catalog #: PEH9270



For Research Use Only , Not for Diagnostic Use