


Product Name: Recombinant Human SCYL3 Protein
Catalog #: PEH9274



Summary

Name	PACE1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8IZE3
Host	E.coli
Species	Human
Molecular Mass	13 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	DADQIEPPKV SSQERPLKVP SELGLGEEFT IQVKKKPVKD PEMDWFADMI PEIKPSAAFL ILPELRTEMV PKKDDVSPVM QFSSKFAAAE ITEGEAEGWE

image



Background

Alternative Names Protein-associating with the carboxyl-terminal domain of ezrin (Ezrin-binding protein PACE-1) (SCY1-like protein 3)

Note

For Research Use Only , Not for Diagnostic Use