


Product Name: Recombinant Human AHI1 Protein
Catalog #: PEH9277



Summary

Name	AHI1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8N157
Host	E.coli
Species	Human
Molecular Mass	35 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	AYQCHVTEEM AKEIKRKIRK KLKEQLTYFP SDTLFHDDKL SSEKRKKKKE VPVFSKAETS TLTISGDTVE GEQKKESSVR SVSSDSHQDD EISSMEQSTE DSMQDDTKPK PKKTKKKTKA VADNNEDVDG DGVHEITSRD SPVYPKCLLD DDLVLGVYIH RTDRLKSDFM ISHPMVKIHV VDEHTGQYVK KDDSGRPVSS YYEKENVDI LPIMTQPYDF KQLKSRLPEW EEQIVFNENF PYLLRGSDES PKVILFFEIL DFLSVDEIKN NSEVQIQECCG

image



Background

Alternative Names Jouberin (Abelson helper integration site 1 protein homolog) (AHI-1)

Product Name: Recombinant Human AHI1 Protein
Catalog #: PEH9277



Note

For Research Use Only , Not for Diagnostic Use