


Product Name: Recombinant Human NDUFAF2 Protein
Catalog #: PEH9278



Summary

Name	NDUF2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8N183
Host	E.coli
Species	Human
Molecular Mass	17 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	YYYYIPQYKN WRGQTIREKR IVEAANKKEV DYEAGDIPTE WEAWIRTRK TPPTMEEILK NEKHREEIKI KSQDFYEKEK LLSKETSEEL LPPPVQTQIK GHASAPYFGK EEPVAPSST GKTFQPGSWM

image



Background

Alternative Names	NADH dehydrogenase [ubiquinone] 1 alpha subcomplex assembly factor 2 (B17.2-like) (B17.2L) (Mimitin) (Myc-induced mitochondrial protein) (MMTN) (NDUFA12-like protein)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human NDUFAF2 Protein
Catalog #: PEH9278

