


Product Name: Recombinant Human NME8 Protein
Catalog #: PEH9281



Summary

Name	TXND3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8N427
Host	E.coli
Species	Human
Molecular Mass	47 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	PEIPLVSDSDS EVSEESPCES VQELYSIAII KPDAVISKKV LEIKRKITKA GFIIIEAEHKT VLTEEQVVNF YSRIADQCDF EEFVSFMTSG LSYILVVSQG SKHNPPSEET EPQTDTEPNE RSEDQPEVEA QVTPGMMKNK QDSLQEYLER QHLAQLCDIE EDAANVAKFM DAFFPDFKKM KSMKLEKTLA LLRPNLFHER KDDVLRRIKD EDFKILEQRQ VVLSEKEAQA LCKEYENEDY FNKLIENMTS GPSLALVLLR DNGLQYWQQL LGPRTVEEAI EYFPESLCAQ FAMDSLVPVNO LYGSDSLETA EREIQHFFPL QSTLGLIKPH ATSEQREQIL KIVKEAGFDL TQVKKMFLTP EQIEKIYPKV TGKDFYKDLL EMLSVGPSMV

image



Background

Alternative Names	Thioredoxin domain-containing protein 3 (NM23-H8) (NME/NM23 family)
--------------------------	---

Product Name: Recombinant Human NME8 Protein
Catalog #: PEH9281



member 8) (Spermatid-specific thioredoxin-2) (Sptrx-2)

Note

For Research Use Only , Not for Diagnostic Use