


Product Name: Recombinant Human RNF10 Protein
Catalog #: PEH9286



Summary

Name	RNF10_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8N5U6
Host	E.coli
Species	Human
Molecular Mass	23 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	TESHQYVVG D TITMQLMKRE KGV LVALPKS KWMNVDPHPIH LGDEQHSQYS KLLASKEQV LHRVVLEEKV ALEQQLAEEK HTPESCFIEA AIQELKTREE ALSGLAGSRR EVTGVVAAL EQLVLMAPLAK ESVFQPRKGV LEYLSAFDEE TTEVCSLDTP SRPLALPLVE EEEAVSEPEP EGLPEACDDL ELADDNLKEG

image



Background

Alternative Names RING finger protein 10

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human RNF10 Protein
Catalog #: PEH9286

