



Product Name: Recombinant Human KCTD17 Protein
Catalog #: PEH9287



Summary

Name	KCD17_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8N5Z5
Host	E.coli
Species	Human
Molecular Mass	27 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EAGEAAPPAG AGGRAAGGWG KVVRLNVGGT VFLTTRQTLC REQKSFLSRL CQGEELQSDR DETGAYLIDR DPTYFGPILN FLRHGKLVLD KDMAEEGVLE EAEFYNIGPL IRIIKDRMEE KDYTVTQVPP KHVYRVLQCQ EEELTQMVST MSDGWRFEQL VNIGSSYNYG SEDQAEFLCV VSKELHSTPN GLSSESSRKT KSTEEQLEEQ QQQEEEEVEEV EVEQVQVEAD

image



Background

Alternative Names BTB/POZ domain-containing protein KCTD17

Note

Product Name: Recombinant Human KCTD17 Protein
Catalog #: PEH9287



For Research Use Only , Not for Diagnostic Use