


Product Name: Recombinant Human C7orf31 Protein
Catalog #: PEH9294



Summary

Name	CG031_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8N865
Host	E.coli
Species	Human
Molecular Mass	24 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KNQDIEEKII EEQSLLSTYE LPSCYPTKDL TSIYDIKFPF KITDTKKTED LYWRQQLSKT QPTYCKPDH WIHYENLKSP LRDQYNMCPD PVSLKPSVL QNKQDTEAFT LEHFLSKPEE ELFLNMENNE ETRPVLGWIP RAGVTKPQTN LLELKNSFSK TGAQKRFHKS ILEDHKDLRD NEHSGMKHQF

image



Background

Alternative Names Uncharacterized protein C7orf31

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human C7orf31 Protein
Catalog #: PEH9294

