


Product Name: Recombinant Human ENDOV Protein
Catalog #: PEH9295



Summary

Name	ENDO_V_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8N8Q3
Host	E.coli
Species	Human
Molecular Mass	25 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EETLSLWKRE QARLKAHVVD RDTEAWQRDP AFSGLQRVGG VDVSFVKGDS VRACASLVVL SFPELEVVE ESRMVSLTAP YVSGFLAFRE VPFLLELVQQ LREKEPGLMP QVLLVDGNGV LHRGFGVAC HLGVLTDLPC VGVAKLLQV DGLENNALHK EKIRLLQTRG DSFPLLGDSG TVLGMALRSH DRSTRPLYIS VGHRMSLEAA

image



Background

Alternative Names Endonuclease V (hEndoV) (EC 3.1.26.-) (Inosine-specific endoribonuclease)

Note

Product Name: Recombinant Human ENDOV Protein
Catalog #: PEH9295



For Research Use Only , Not for Diagnostic Use