


Product Name: Recombinant Human CFAP43 Protein
Catalog #: PEH9301



Summary

Name	CFA43_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8NDM7
Host	E.coli
Species	Human
Molecular Mass	31 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KENHYLSTTP KVSVDLGSDS EHTKQKASTD LSQDELVLTD VKKEIPWIIQQ KSQEAIKKEV NLFSSKKRKEI KQGIKSLSKT ILNMMEENDK LENIAKLDQQ EFGDLLEELE RLHDESQEEV AKMIKDVEMH NLAKSYLAEL IKEECWNSMA VKGRALKCFH IPCVVNFPM KARTVEELKE LERVLQKKI EAELKLRKE IVEAQSGVKL IKQRHEEDDE EEEEEDKTVK YSNLPNYLLG SLSTDFGVDT SLLSSQLELH

image



Background

Alternative Names	Cilia- and flagella-associated protein 43 (WD repeat-containing protein 96)
--------------------------	---

Product Name: Recombinant Human CFAP43 Protein
Catalog #: PEH9301



Note

For Research Use Only , Not for Diagnostic Use