


Product Name: Recombinant Human DUSP18 Protein
Catalog #: PEH9303



Summary

Name	DUS18_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8NEJ0
Host	E.coli
Species	Human
Molecular Mass	20 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SQITKSLYIS NGVAANNKLM LSSNQITMVI NVSVEVVNTL YEDIQYMQVP VADSPNSRLC DFFDPIADHI HSVEMKQGRT LLHCAAGVSR SAALCLAYLM KYHAMSL LDA HTWTKSCRPI IRPNSGFWEQ LIHYEFQLFG KNTVHVMVSSP VGMIPDIYEK

image



Background

Alternative Names	Dual specificity protein phosphatase 18 (EC 3.1.3.16) (EC 3.1.3.48) (Low molecular weight dual specificity phosphatase 20) (LMW-DSP20)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human DUSP18 Protein
Catalog #: PEH9303

