



Product Name: Recombinant Human AGBL3 Protein
Catalog #: PEH9305



Summary

Name	CBPC3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8NEM8
Host	E.coli
Species	Human
Molecular Mass	12 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs.Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	QSKPHRSLES LSPLKGPKKN KHSQIWAIGN EDIKPLSSKW ETASSSFGMD ANVLKYKSLQ AEETNQQSSK HTALHLTKNK DEQANKNDGQ

image



Background

Alternative Names Cytosolic carboxypeptidase 3 (EC 3.4.17.-) (ATP/GTP-binding protein-like 3)

Note

For Research Use Only , Not for Diagnostic Use