


Product Name: Recombinant Human NAV1 Protein
Catalog #: PEH9306



Summary

Name	NAV1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8NEY1
Host	E.coli
Species	Human
Molecular Mass	12 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	DELLSSKAKA QKSSGPVPSA KGQEERAFK VDPELVVTVL GDLEQLLFSQ MLDPESQRKR TVQNVLDLRQ NLEETMSSLR GSQVTHSSLE

image



Background

Alternative Names Neuron navigator 1 (Pore membrane and/or filament-interacting-like protein 3) (Steerin-1) (Unc-53 homolog 1) (unc53H1)

Note

For Research Use Only , Not for Diagnostic Use