

Product Name: Recombinant Human ST13P5 Protein
Catalog #: PEH9307



Summary

Name	F10A5_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8NFI4
Host	E.coli
Species	Human
Molecular Mass	42 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	AFVKMCKQDP SVLHTEEMRF LREWVESMGG KVPPATQKAK SEENTKEEKP DSKKVEEDLK ADEPSTEESD LEIDKEGVIE PDTDAPQEMG DENVEITEEM MDQANDKKVA AIEVLNDGEL QKAIDLFTDA IKLNPRLAIL YAKRASVFK LQKPNAAIQD CDRAIEINPD SAQPYKWRGK AHRLLGHWEE AAHDLAFACK LDYDEDASAM LKEVQPRAQK IAEHWRKYER KHEEREIKER IERVKKAQEE QERAQREEEA RRQSGAHYGP FPGGFPGGMP GNFPGGMPGM GGDMPGMAGM PGLNEILSDP EALAAMQDPE VMVAFQDVAQ NPANMSKYQS NPKVMNLISK

image



Background

Alternative Names	Putative protein FAM10A5 (Suppression of tumorigenicity 13 pseudogene 5)
--------------------------	--

Product Name: Recombinant Human ST13P5 Protein
Catalog #: PEH9307



Note

For Research Use Only , Not for Diagnostic Use