

Product Name: Recombinant Human BBS1 Protein
Catalog #: PEH9308



Summary

Name	BBS1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8NFJ9
Host	E.coli
Species	Human
Molecular Mass	13 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	QAKEDRIDPL TLKEMLESIR ETAEELSIQ SLRFLQLELS EMEAFVNQHK SNSIKRQTVI TTMTTLKKNL ADEDAVSCLV LGTENKELLV

image



Background

Alternative Names Bardet-Biedl syndrome 1 protein (BBS2-like protein 2)

Note

For Research Use Only , Not for Diagnostic Use