

Product Name: Recombinant Human TNIP2 Protein
Catalog #: PEH9310



Summary

Name	TNIP2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8NFZ5
Host	E.coli
Species	Human
Molecular Mass	11 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KINDCAEVKQ ELAASRTARD AALERVQMLE QQILAYKDDF MSERADRERA QSRIQELEEK VASLLHQVSW RQDSREPDAG

image



Background

Alternative Names	TNFAIP3-interacting protein 2 (A20-binding inhibitor of NF-kappa-B activation 2) (ABIN-2) (Fetal liver LKB1-interacting protein)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use