


Product Name: Recombinant Human GIMAP7 Protein
Catalog #: PEH9312



Summary

Name	GIMA7_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8NHV1
Host	E.coli
Species	Human
Molecular Mass	33 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	IVLVGKTGSG KSATANTILG EEIFDSRIAA QAVTKNCQKA SREWQGRDLL VVDTPGLFDT KESLDTTCKE ISRCIISCP GPHAIVLVLL LGRYTEEEQK TVALIKAVFG KSAMKHMVIL FTRKEELEGQ SFHDFIADAD VGLKSIVKEC GNRCCAFSNS KKTSKAEKES QVQELVELIE KMQVCNEGAY FSDDIYKDTE ERLKQREEVL RKIYTDQLNE EIKLVEEDKH KSEEEKEKEI KLLKLKYDEK IKNIREEAER NIFKDVFNRI

image



Background

Alternative Names	GTPase IMAP family member 7 (Immunity-associated nucleotide 7 protein) (IAN-7)
--------------------------	--

Product Name: Recombinant Human GIMAP7 Protein
Catalog #: PEH9312



Note

For Research Use Only , Not for Diagnostic Use