


Product Name: Recombinant Human PI4K2B Protein
Catalog #: PEH9315



Summary

Name	P4K2B_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8TCG2
Host	E.coli
Species	Human
Molecular Mass	13 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	ASADGGSP EEEDGEREPLL PRIAWAHPRR GAPGSAVRLL DAAGEEGEAG DEELPLPPGD VGVSRSSAE LDRSRPAVSV TIGTSEMNAF LDDPEFADIM

image



Background

Alternative Names	Phosphatidylinositol 4-kinase type 2-beta (EC 2.7.1.67) (Phosphatidylinositol 4-kinase type II-beta) (PI4KII-BETA)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use