


Product Name: Recombinant Human CABLES1 Protein
Catalog #: PEH9317



Summary

Name	CABL1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8TDN4
Host	E.coli
Species	Human
Molecular Mass	32 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	FSVLPYRDST QVGD LKLDGG RQSTGAVSLK EIIGLEGVEL GADGKTVSYT QFLPTNAFG ARRNTIDSTS SFSQFRNL SH RSLSIGRASG TQGS LDTGSD LGDFMDYDPN LLDDPQWPCG KHKRVLI FPS YMTTVIDYVK PSDLKKDMNE TFKEKFP HIK LTL SKIRSLK REMRKLAQED CGLE EPTVAM AFVYFEKLAL KGKLNKQNRK LCAGACVLLA AKIGSDLKKH EVKHLIDKLE EKFR LNRREL IAFEFPVLVA LEFALHLPEH

image



Background

Alternative Names CDK5 and ABL1 enzyme substrate 1 (Interactor with CDK3 1) (Ik3-1)

Product Name: Recombinant Human CABLES1 Protein
Catalog #: PEH9317



Note

For Research Use Only , Not for Diagnostic Use