


Product Name: Recombinant Human OTUD7A Protein
Catalog #: PEH9319




Summary

Name	OTU7A_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8TE49
Host	E.coli
Species	Human
Molecular Mass	25 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	APLAQPESPT ASAGEDVQSL ADSLDSRDS VCSNSNSNNG KNGKDKEKEK QRKEKDKTRA DSVANKLGSF SKTLGIKLLK NMGGLGGLVH GKMGRANSAN GKNGDSAERG KEKKAKSRKG SKEESGASAS TSPSEKTTPS PTDKAAGASP AEKGGGPRGD AWKYSTDVKL SLNILRAAMQ GERKFIFAGL LLTSHRHQFH EEMIGYYLTS

image



SDS



WB

Background

Alternative Names	OTU domain-containing protein 7A (EC 3.4.19.12) (Zinc finger protein Cezanne 2)
--------------------------	---

Note

Product Name: Recombinant Human OTUD7A Protein
Catalog #: PEH9319



For Research Use Only , Not for Diagnostic Use