



Product Name: Recombinant Human PSPC1 Protein
Catalog #: PEH9326



Summary

Name	PSPC1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8WXF1
Host	E.coli
Species	Human
Molecular Mass	44 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MGFTIDIKSF LKPGEKTYTQ RCRLFVGNLP TDITEEDFKR LFERYGEPSE VFINRDRGFG FIRLESRTLA EIAKAELDGT ILKSRLRIR FATHGAALTV KNLSPVVSNE LLEQAFSQFG PVEKAVVVVD DRGRATGKGF VEFAAKPPAR KALERCGDGA FLLTTTPRPV IVEPMEQFDD EDGLPEKLMQ KTQQYHKERE QPPRFAQPGT FEFYASRWK ALDEMEKQQR EQVDRNIREA KEKLEAEMEA ARHEHQLMLM RQDLMRRQEE LRRLEELRNQ ELQKRKQIQL RHEEEHRRRE EEMIRHREQE ELRRQQEGFK PNYMENREQE MRMGDMGPRG AINMGDAFSP

image



Background

Alternative Names Paraspeckle component 1 (Paraspeckle protein 1)

Product Name: Recombinant Human PSPC1 Protein
Catalog #: PEH9326



Note

For Research Use Only , Not for Diagnostic Use