


Product Name: Recombinant Human CTNNBL1 Protein
Catalog #: PEH9328



Summary

Name	CTBL1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8WYA6
Host	E.coli
Species	Human
Molecular Mass	18 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	PNRGTKRPRD DEEEQKMRR KQTGTRERGR YREEEMTVVE EADDDKKRLL QIIDRDGEEE EEEEEPLDES SVKKMILTFE KRSYKNQELR IKFPDNPEKF MESELDLNDI IQEMHVVATM PDLYHLLVEL NAVQSLLGLL

image



Background

Alternative Names	Beta-catenin-like protein 1 (Nuclear-associated protein) (NAP) (Testis development protein NYD-SP19)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use