


Product Name: Recombinant Human ACOT12 Protein
Catalog #: PEH9329



Summary

Name	ACO12_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8WYK0
Host	E.coli
Species	Human
Molecular Mass	20 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MSQAIQPAHA TARGELSAGQ LLKWIDTTAC LAAEKHAGVS CVTASVDDIQ FEETARVGQV ITIKAKVTRA FSTSMEISIK VMVQDMLTGI EKLVSVAFST FVAKPVGKEK IHLKPVTLLT EQDHVEHNLA AERRKVRLQH EDTFNNLMKE SSKFDDLIFD

image



Background

Alternative Names	Acetyl-coenzyme A thioesterase (EC 3.1.2.1) (Acyl-CoA thioester hydrolase 12) (Acyl-coenzyme A thioesterase 12) (Acyl-CoA thioesterase 12) (Cytoplasmic acetyl-CoA hydrolase 1) (CACH-1) (hCACH-1) (START domain-containing protein 15) (StARD15)
--------------------------	---

Product Name: Recombinant Human ACOT12 Protein
Catalog #: PEH9329



Note

For Research Use Only , Not for Diagnostic Use