



Product Name: Recombinant Human AHCTF1 Protein
Catalog #: PEH9330



Summary

Name	ELYS_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8WYP5
Host	E.coli
Species	Human
Molecular Mass	14 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SQNIMCPTKL VKSAFKTAQE TSTMTMNVSQ VDDVVSSKTR TRGQRIQNVN VKSAQQEASA DVATPKMPGQ SVRKKTRKAK EISEASENIY SDVRGLSQNQ QIPQNSVTPR

image



Background

Alternative Names	Protein ELYS (Embryonic large molecule derived from yolk sac) (Protein MEL-28) (Putative AT-hook-containing transcription factor 1)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use