



Product Name: Recombinant Human RAPGEF4 Protein
Catalog #: PEH9332



Summary

Name	RPGF4_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8WZA2
Host	E.coli
Species	Human
Molecular Mass	23 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	TEEEKKECDE ELQDTMLLS QMGPDAMMRM ILRKPPGQRT VDDLEIIYEE LLHIKALSHL STTVKRELAG VLIFESHAKG GTVLFNQGEE GTSWYIILKG SVNVVIYGKG VVCTLHEGDD FGKLALVND A PRAASIVLRE DNCHFLRVDK EDFNRLRDV EANTVRLKEH DQDVLVLEKV PAGNRASNQG

image



Background

Alternative Names	Rap guanine nucleotide exchange factor 4 (Exchange factor directly activated by cAMP 2) (Exchange protein directly activated by cAMP 2) (EPAC 2) (cAMP-regulated guanine nucleotide exchange factor II) (cAMP-GEFII)
--------------------------	--

Note

Product Name: Recombinant Human RAPGEF4 Protein
Catalog #: PEH9332



For Research Use Only , Not for Diagnostic Use