


Product Name: Recombinant Human DLG3 Protein
Catalog #: PEH9335



Summary

Name	DLG3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q92796
Host	E.coli
Species	Human
Molecular Mass	59 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SGPSWWPECT CTNRDWYEQV NGSDGMFKYE EIVLERGNSG LGFSIAGGID NPHVPDDPGI FITKIIPGGA AAMDGRLGVN DCVLRVNEVD VSEVVHSRAV EALKEAGPVV RLVVRRRQPP PETIMEVNLL KGPKGLGFSI AGGIGNQHIP GDNSIYITKI IEGGAAQKDG RLQIGDRLLA VNNTNLQDVR HEEAVASLKN TSDMVYLKVA KPGSLHLNDM YAPPDYASTF TALADNHISH NSSLGYLGA V ESKVSYPAPP QVPTRYSPI PRHMLAEEDF TREPRKIILH KGSTGLGFNI VGGEDGEGIF VSFILAGGPA DLSGELRRGD RILSVNGVNL RNATHEQAAA ALKRAGQSVT IVAQYRPEEY SRFESKIHDL REQMMNSSMS SGSGSLRTSE KRSLYVRALF DYDRTRDSC L PSQGLSFSYG DILHVINASD DEWWQARLVT PHGESEQIGV IPSKKRVEKK ERARLKTVKF HARTGMIESN RDFPGLSDDY YGAKNLKGQE DAILYSEPV T

image



Background

Product Name: Recombinant Human DLG3 Protein
Catalog #: PEH9335



Alternative Names

Disks large homolog 3 (Neuroendocrine-DLG) (Synapse-associated protein 102) (SAP-102) (SAP102) (XLMR)

Note

For Research Use Only , Not for Diagnostic Use