


Product Name: Recombinant Human PMM1 Protein
Catalog #: PEH9338



Summary

Name	PMM1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q92871
Host	E.coli
Species	Human
Molecular Mass	31 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MAVTAQAARR KERVLCFLDV DGLTPARQK IDPEVAAFQ KLSRVQIGV VGGSDYCKIA EQLGDGDEVI EKFDYVFAEN GTVQYKHGRL LSKQTIQNH GEELLQDLIN FCLSYMALLR LPKKGTFIE FRNGMLNISP IGRSCTLEER IEFSELDKKE KIREKFVEAL KTEFAGKGLR FSRGGMISFD VFPEGWDKRY CLDSLQDSF DTIHFFGNET SPGGNDFEIF ADPRTVGHSV VSPQDTVQRC REIFFPETAH EA

image



Background

Alternative Names Phosphomannomutase 1 (PMM 1) (EC 5.4.2.8) (PMMH-22)

Product Name: Recombinant Human PMM1 Protein
Catalog #: PEH9338



Note

For Research Use Only , Not for Diagnostic Use