



Summary

Name	ARHG2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q92974
Host	E.coli
Species	Human
Molecular Mass	24 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	DRALVELLRE KVGLEAEMTH FQAEEDGGSG MALPTLPRGL FRSESLESPR GERLLQDAIR EVEGLKDLLV GPGVELLLTP REPALPLEPD SGGNTSPGVT ANGEARTFNG SIELCRADSD SSQRDRNGNQ LRSPQEEALQ RLVNLYGLLH GLQAAVAQQD TLMEARFPEG PERREKLCRA NSRDGEAGRA GAAPVAPEKQ

image



Background

Alternative Names	Rho guanine nucleotide exchange factor 2 (Guanine nucleotide exchange factor H1) (GEF-H1) (Microtubule-regulated Rho-GEF) (Proliferating cell nucleolar antigen p40)
--------------------------	--

Note

Product Name: Recombinant Human ARHGEF2 Protein
Catalog #: PEH9339



For Research Use Only , Not for Diagnostic Use