

Product Name: Recombinant Human UBE2E3 Protein
Catalog #: PEH9343



Summary

Name	UB2E3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q969T4
Host	E.coli
Species	Human
Molecular Mass	10 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	ESPSTSSGSS DADQRDPAAP EPEEQEERKP SATQQKKNTK LSSKTTAKLS TSAKRIQKEL AEITLDPPPN CSAGPKGDNI

image



Background

Alternative Names	Ubiquitin-conjugating enzyme E2 E3 (EC 2.3.2.23) (E2 ubiquitin-conjugating enzyme E3) (UbcH9) (Ubiquitin carrier protein E3) (Ubiquitin-conjugating enzyme E2-23 kDa) (Ubiquitin-protein ligase E3)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use