

Product Name: Recombinant Human ATG4C Protein
Catalog #: PEH9351



Summary

Name	ATG4C_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q96DT6
Host	E.coli
Species	Human
Molecular Mass	30 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	AESGCTIEDH VIAGNVEEFR KDFISRIWLT YREEFPQIEG SALTDCGWG CTLRTGQMLL AQGLILHFLG RAWTWPDALN IENSSESWT SHTVKKFTAS FEASLSGERE FKTPTISLKE TIGKYSDDE MRNEVYHRKI ISWFGDSPLA LFGLHQLIEY GKKS GKKAGD WYGPVVAHI LRKAVEEARH PDLQGITYV AQDCTVYNSD VIDKQSASMT SDNADDKAVI ILVPVRLGGE

image



Background

Alternative Names	Cysteine protease ATG4C (EC 3.4.22.-) (AUT-like 3 cysteine endopeptidase) (Autophagin-3) (Autophagy-related cysteine endopeptidase 3) (Autophagy-related protein 4 homolog C)
--------------------------	---

Product Name: Recombinant Human ATG4C Protein
Catalog #: PEH9351



Note

For Research Use Only , Not for Diagnostic Use