

Product Name: Recombinant Human NAF1 Protein
Catalog #: PEH9357



Summary

Name	NAF1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q96HR8
Host	E.coli
Species	Human
Molecular Mass	28 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SSSSCISLP PVLSDGDDDL QIEKENKNFP LTKDELLN ELPSVEELTI ILPEDIELKP LGMVSSIEQ LVIESMTNL PPVNEETVIF KSDRQAAGKI FEIFGPVAHP FYVLRFNSSD HIESKGIKIK ETMYFAPSMK DFTQYIFTEK LKQDKGSDAS WKNDQEPPPE ALDFSDDKE KEAKQRKKSQ IQGRKKLKSE FNEPGEDFTE VHQNWNAHSS ASEHAKGYRN

image



Background

Alternative Names H/ACA ribonucleoprotein complex non-core subunit NAF1 (hNAF1)

Note

Product Name: Recombinant Human NAF1 Protein
Catalog #: PEH9357



For Research Use Only , Not for Diagnostic Use