


Product Name: Recombinant Human ZFYVE19 Protein
Catalog #: PEH9360



Summary

| | |
|------------------------------|--|
| Name | ANCHR_HUMAN |
| Purity | Greater than 85% as determined by reducing SDS-PAGE |
| Endotoxin level | Please contact with the lab for this information |
| Tag | C-His |
| Accession # | Q96K21 |
| Host | E.coli |
| Species | Human |
| Molecular Mass | 11 kDa |
| Buffer | Phosphate buffered saline |
| Shipping | The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below. |
| Stability&Storage | Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles. |
| Sequence | TRQDQMIAER LARLRQENKP KLVPSQAEIE ARLAALKDER QGSIPSTQEM EARLAALQGR VLPSQTPQPA HHTPDTRTQA |

image



Background

| | |
|--------------------------|--|
| Alternative Names | Abscission/NoCut checkpoint regulator (ANCHR) (MLL partner containing FYVE domain) (Zinc finger FYVE domain-containing protein 19) |
|--------------------------|--|

Note

For Research Use Only , Not for Diagnostic Use