


Product Name: Recombinant Human PHACTR3 Protein
Catalog #: PEH9362



Summary

Name	PHAR3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q96KR7
Host	E.coli
Species	Human
Molecular Mass	15 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SSVERGKERE EAWSFDGALE NKRTAAKESE ENKENLIINS ELKDDLILLYQ DEEALNDSII SGTLPKCKK ELLAVKLNRN PSKQELEDNRN IFPRRTDEER QEIRQQIEMK

image



Background

Alternative Names	Phosphatase and actin regulator 3 (Scaffold-associated PP1-inhibiting protein) (Scapinin)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use