


## Summary

<b>Name</b>	NALP4_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	C-His
<b>Accession #</b>	Q96MN2
<b>Host</b>	E.coli
<b>Species</b>	Human
<b>Molecular Mass</b>	16 kDa
<b>Buffer</b>	Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	ISNQNLKILQ IGCNEIGDVG VQLLCRALTH TDCRLEILGL EECGLTSTCC KDLASVLTCS KTLQQNLNLT NTLDHTGVVV LCEALRHPEC ALQVLGLRKT DFDEETQALL TAEERNPNL

## image



## Background

<b>Alternative Names</b>	NACHT, LRR and PYD domains-containing protein 4 (Cancer/testis antigen 58) (CT58) (PAAD and NACHT-containing protein 2) (PAN2) (PYRIN and NACHT-containing protein 2) (PYRIN-containing APAF1-like protein 4) (PYPAF4) (Ribonuclease inhibitor 2)
--------------------------	---

## Note

**Product Name: Recombinant Human NLRP4 Protein**  
**Catalog #: PEH9367**



---

For Research Use Only , Not for Diagnostic Use