


Product Name: Recombinant Human JAKMIP1 Protein
Catalog #: PEH9369



Summary

Name	JKIP1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q96N16
Host	E.coli
Species	Human
Molecular Mass	50 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KPEMETDAVQ MANEELRAKL TSIQIEFQQE KSKVGKLRER LQEAKLEREQ EQRRHTAYIS ELKAKLHEEK TKELQALREG LIRQHEQEAA RTAKIKEGEL QRLQATLNVL RDGAADKVKT ALLTEAREEA RRAFDGERLR LQQEILELKA ARKQAEELALS NCMQADKTKA ADLRAAYQAH QDEVHRIKRE CERDIRRLMD EIKGKDRVIL ALEKELGVQA GQTQKLLQK EALDEQLVQV KEAERHHSSP KRELPPGIGD MVELMGVQDQ HMDERDVRRF QLKIAELNSV IRKLEDNRTL LADERNELLK RSRETEVQLK PLVEKNKRMN KKNEDLLQSI QRMEEKIKNL TRENDEMKEK LSAQASLKRH TSLNDLSLTR DEQEIEFLRL QVLEQQHVID DLSLERERLL

image



Background

Product Name: Recombinant Human JAKMIP1 Protein
Catalog #: PEH9369



Alternative Names

Janus kinase and microtubule-interacting protein 1 (GABA-B receptor-binding protein) (Multiple alpha-helices and RNA-linker protein 1) (Marlin-1)

Note

For Research Use Only , Not for Diagnostic Use